

EYE SURGERY



CATALOGUE



Eye

1- Scissors.....	01-17
2- Micro Needle Holders.....	18-20
3- Forceps.....	21-33
4- Specula.....	34-36
5- Hooks and Retractors.....	37-38
6- Hooks / Retractors / Manipulators.....	39-50
7- Fixation Rings and Markers.....	51-53
8- Probes and Dilators.....	54-55
9- Lacrimal Intubation Sets.....	56-57
10- Cannulas.....	58-64
11- Miscellaneous.....	65


Eye Scissors




EL-4767-EL-4772




EL-4767 $\frac{1}{1}$




EL-4769 $\frac{1}{1}$



EL-4768 $\frac{1}{1}$




EL-4773-EL-4774




EL-4773 $\frac{1}{1}$




EL-4770 $\frac{1}{1}$



EL-4771 $\frac{1}{1}$







EL-4772 $\frac{1}{1}$



EL-4774 $\frac{1}{1}$

Iris Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4767	S/S	STR	RIBBON	4
EL-4768	S/S	STR	RIBBON	4 1/4
EL-4769	B/B	STR	RIBBON	4 1/4
EL-4770	S/S	CVD	RIBBON	4
EL-4771	S/S	CVD	RIBBON	4 1/4
EL-4772	B/B	CVD	RIBBON	4 1/4
EL-4773	S/S	STR	SQUARE	4 1/2
EL-4774	S/S	CVD	SQUARE	4 1/2





Iris Scissors (continued)


Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4775	S/S	STR	ROUND	3 1/2
EL-4776	S/S	STR	ROUND	3 3/4
EL-4777	S/S	STR	ROUND	4 1/8
EL-4778	S/S	STR	ROUND	4 1/4
EL-4779	S/B	STR	ROUND	4 1/4
EL-4780	S/S	STR	ROUND	4 1/2
EL-4781	S/S	STR	ROUND	4 3/4
EL-4782	S/S	CVD	ROUND	3 1/2
EL-4783	S/S	CVD	ROUND	3 3/4
EL-4784	S/S	CVD	ROUND	4 1/8
EL-4785	S/S	CVD	ROUND	4 1/4
EL-4786	S/B	CVD	ROUND	4 1/4
EL-4787	S/S	CVD	ROUND	4 1/2
EL-4788	S/S	CVD	ROUND	4 3/4
EL-4789	S/S	LAT-ANG	ROUND	4 1/4
EL-4790	S/S	STR	ROUND, CENTERED	4 1/4




Eye Scissors



EL-4792-EL-4794



EL-4792 1/1



EL-4793 1/1




EL-4794 1/1




EL-4791-EL-4792


1/2



EL-4791 1/1




EL-4795 1/1




EL-4796

1/2




EL-4796 1/1




EL-4797-EL-4800


1/2




EL-4797 1/1



EL-4799 1/1




EL-4798 1/1



EL-4800 1/1


STEVENS Tenotomy Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4791	S/S	STR	RIBBON	4
EL-4792	B/B	STR	RIBBON	3 3/4
EL-4793	B/B	STR	RIBBON	3 7/8
EL-4794	B/B	CVD	RIBBON	3 3/4
EL-4795	S/S	CVD	RIBBON	4
EL-4796	B/B	CVD	RIBBON	4 1/4
EL-4797	S/S	STR	SQUARE	4 1/2
EL-4798	B/B	STR	SQUARE	4 1/2
EL-4799	S/S	CVD	SQUARE	4 1/2
EL-4800	B/B	CVD	SQUARE	4 1/2




EL-4801-EL-4806

1/2




EL-4801

1/1




EL-4805

1/1




EL-4802

1/1




EL-4806

1/1




EL-4803-EL-4807

1/2




EL-4803

1/1




EL-4807

1/1




EL-4804-EL-4808

1/2



EL-4804

1/1



EL-4808

1/1



EL-4809

1/2



EL-4809

1/1

**STEVENS Tenotomy Scissors
(continued)**












Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4801	S/S	STR	ROUND	4 1/4
EL-4802	B/B	STR	ROUND	4 1/4
EL-4803	B/B	STR	ROUND	5
EL-4804	B/B	STR	ROUND	5 3/4
EL-4805	S/S	CVD	ROUND	4 1/4
EL-4806	B/B	CVD	ROUND	4 1/4
EL-4807	B/B	CVD	ROUND	5
EL-4808	B/B	CVD	ROUND	5 3/4

JAMESON-WERBER Tenotomy Scissor

Catalog Number	Tip	Blade	Length (in.)
EL-4809	B/B	DEL, CVD	5 1/8



Eye Scissors



Strabismus Scissors, Delicate

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4810	B/B	CVD	RIBBON	4 1/4
EL-4811	B/B	STR	SQUARE	4 1/2
EL-4812	B/B	CVD	SQUARE	4 1/2
EL-4813	B/B	STR	ROUND	4 1/2
EL-4814	B/B	CVD	ROUND	4 1/2

KNAPP Strabismus Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
EL-4815	B/B	STR	4 3/4
EL-4816	B/B	CVD	4 3/4



EL-4817-EL-4818

Enucleation Scissors

Catalog Number	Tip	Blade	Length (in.)
EL-4817	B/B	SLIGHT CVD	5
EL-4818	B/B	MEDIUM CVD	5 1/4
EL-4819	B/B	FULL CVD	5 1/8



EL-4817



EL-4818




EL-4819




EL-4820

Eye Scissors

Catalog Number	Tip	Blade	Length (in.)
EL-4820	PROBE	ANG	4 1/2



EL-4820




EL-4821-EL-4823

Utility Scissors, Ribbon Handle


Catalog Number	Tip	Blade	Length (in.)
EL-4821	B/B	STR	4
EL-4822	B/B	CVD	4
EL-4823	B/B	CVD	4



EL-4821



EL-4822




EL-4823




Eye Scissors




EL-4824-EL-4825



EL-4824




EL-4825



EL-4826



EL-4826



EL-4827



EL-4827



EL-4828



EL-4828

Stitch Scissors, Ribbon Handle

Catalog Number	Tip	Blade	Length (in.)
EL-4824	S/S	CVD	3 3/4
EL-4825	S/S	CVD	4 3/4

Stitch Scissors

Catalog Number	Tip	Blade	Length (in.)
EL-4826	S/S	STR	4

WILMER Conjunctival Scissors


Catalog Number	Tip	Blade	Length (in.)
EL-4827	S/S	ANG-FLAT	4 1/8

IOL Cutter


Catalog Number	Blade	Blade Length (mm)	Diameter (mm)	Length (in.)
EL-4828	SERR	6	2	4 1/2

Tungsten Carbide (TC) Scissors


- Excellent cutting accuracy
- Long working life
- Gold-plated rings



EL-4829-EL-4830




EL-4829




EL-4830

TC Iris Scissors, Delicate


Catalog Number	Tip	Blade	Length (in.)
EL-4829	S/S	STR	4 1/4
EL-4830	S/S	CVD	4 1/4



EL-4831-EL-4832



EL-4831



EL-4832

TC Strabismus Scissors, Delicate

Catalog Number	Tip	Blade	Length (in.)
EL-4831	B/B	STR	4 1/2
EL-4832	B/B	CVD	4 1/2




Eye TC Scissors





Wavecut

Wavecut scissors have one serrated cutting edge to help prevent tissue from slipping. The **W** after the item number indicates the scissors are the wavecut style, e.g. .




TC Iris Scissors, Wavecut

Catalog Number	Tip	Blade	Length (in.)
EL-4833	S/S	CVD	4 1/4



TC Strabismus Wavecut Scissors

Catalog Number	Tip	Blade	Length (in.)
EL-4834	B/B	CVD	4 1/2





EL-4835-EL-4837



EL-4835



EL-4836



EL-4837

VANNAS Micro Scissors


Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4835	S/S	STR	FLAT, SERR	3 1/4
EL-4836	S/S	CVD	FLAT, SERR	3 1/4
EL-4837	S/S	ANG	FLAT, SERR	3 1/4




EL-4838-EL-4840



EL-4838



EL-4839




EL-4840

VANNAS Micro Scissors


Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4838	S/S	STR	FLAT, SERR	4 3/4
EL-4839	S/S	CVD	FLAT, SERR	4 3/4
EL-4840	S/S	ANG	FLAT, SERR	4 3/4
EL-4841	S/S	STR	FLAT, SERR	6 1/4
EL-4842	S/S	CVD	FLAT, SERR	6 1/4



EL-4841-EL-4842




EL-4841



EL-4842



EL-4843













EL-4843

VANNAS Micro Scissors, Wide Handle

Catalog Number	Tip	Tip Size (mm)	Blade	Handle	Length (in.)
EL-4843	S/S	10	CVD	FLAT, WIDE SERR	4 3/4



Eye Scissors








VANNAS Capsulotomy Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4844	S/S	STR	FLAT, SERR	3 1/4
EL-4845	S/S	CVD	FLAT, SERR	3 1/4
EL-4846	S/S	ANG-DOWN	FLAT, SERR	3 1/4
EL-4847	S/S	ANG-UP	FLAT, SERR	3 1/4

Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4848	S/S	STR	FLAT, SERR	4 1/8
EL-4849	S/S	CVD	FLAT, SERR	4 1/8
EL-4850	S/S	STR	FLAT, SERR	4 3/4
EL-4851	S/S	CVD	FLAT, SERR	4 3/4
EL-4852	S/S	STR	FLAT, SERR	6 1/4
EL-4853	S/S	CVD	FLAT, SERR	6 1/4
EL-4854	S/S	STR	ROUND	5 3/4
EL-4855	S/S	CVD	ROUND	5 3/4
EL-4856	B/B	CVD	ROUND	5 3/4



MCCLURE Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4857	S/S	STR	FLAT, SERR	3 3/4

NOYES Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4858	S/B	STR	FLAT, SERR	5

CASTROVIEJO Corneal Scissors


Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4859	B/B	LAT-ANG	FLAT SERR	3 3/4
EL-4860	B/B	LAT-CVD	FLAT SERR	3 3/4
EL-4861	B/B	LEFT-CVD	FLAT SERR	4 3/8
EL-4862	B/B	RIGHT-CVD	FLAT SERR	4 3/8

CASTROVIEJO Corneal Section Scissors


Catalog Number	Tip	Tip Size (mm)	Blade	Handle	Length (in.)
EL-4863	B/B	9	CVD, ANG RIGHT	FLAT, SERR	4 1/4
EL-4864	B/B	9	CVD, ANG LEFT	FLAT, SERR	4 1/4



Eye Scissors



EL-4865



EL-4865




EL-4866



EL-4866



EL-4867



EL-4867



EL-4868



EL-4868



EL-4869



EL-4869



EL-4870



EL-4870



EL-4871



EL-4871

CASTROVIEJO Universal Corneal Scissors

Catalog Number	Tip	Tip Size	Blade	Handle	Length (in.)
EL-4865	B/B	16 MM FROM MID-SCREW	CVD	FLAT, FENESTRATED	4 3/4

CASTROVIEJO Corneal Scissors






Catalog Number	Tip	Tip Size	Blade	Handle	Length (in.)
EL-4866	B/B	SMALL	CVD RIGHT	FLAT, WIDE SERR	4
EL-4867	B/B	SMALL	CVD LEFT	FLAT, WIDE SERR	4

TROUTMAN-CASTROVIEJO Corneal Scissors, W/ Catch

Catalog Number	Tip	Tip Size	Blade	Handle	Length (in.)
EL-4868	B/B	MEDIUM	CVD RIGHT	FLAT, WIDE SERR	4 1/8
EL-4869	B/B	MEDIUM	CVD LEFT	FLAT, WIDE SERR	4 1/8

TROUTMAN-KATZIN Corneal Scissors, W/ Catch

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4870	B/B	STRONG CVD RIGHT	FLAT, WIDE SERR	4 1/8
EL-4871	B/B	STRONG CVD LEFT	FLAT, WIDE SERR	4 1/8



TROUTMAN Corneal Scissors, W/ Catch

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4872	B/B	CVD RIGHT	FLAT, WIDE SERR	4 1/8

KATZIN Corneal Transplant Scissors

Catalog Number	Tip	Tip Size	Blade	Handle	Length (in.)
EL-4873	B/B	MEDIUM	CVD RIGHT	FLAT, FENESTRATED	4 1/8
EL-4874	B/B	MEDIUM	CVD LEFT	FLAT, FENESTRATED	4 1/8

RICH Corneal Scissors





Catalog Number	Tip	Tip Size	Blade	Handle	Length (in.)
EL-4875	S/S	X-DEL	CVD RIGHT	FLAT, WIDE SERR	4 1/8
EL-4876	S/S	X-DEL	CVD LEFT	FLAT, WIDE SERR	4 1/8

Micro Iris Scissors

Catalog Number	Tip	Tip Size	Blade	Handle	Length (in.)
EL-4877	S/S	10 MM FROM MID-SCREW	CVD	FLAT, SERR	4 1/8



Eye Scissors



BARRAQUER Iris Scissors

Catalog Number	Tip	Tip Size (mm)	Blade	Length (in.)
EL-4878	B/B	8	ANG	2 1/4

JAFFE Stitch Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4879	X-SHARP/ X-SHARP	CVD	FLAT, WIDE SERR	4 3/4

Stitch Scissors, Fine

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4880	S/S	STR	FLAT, WIDE SERR	5 1/8

STERN-GILLS Scissors


Catalog Number	Tip	Tip Size (mm)	Blade	Handle	Length (in.)
EL-4881	S/S	10	STR-ANG FORWARD	FLAT, WIDE SERR	4 1/2

URIBE-STERN Scissors

Catalog Number	Tip	Tip Size (mm)	Blade	Handle	Length (in.)
EL-4882	S/S	10	CVD-ANG FORWARD	FLAT, WIDE SERR	4 3/4




EL-4883



EL-4883



EL-4884



EL-4884



EL-4885



EL-4885

UTHOFF-GILLS Capsulotomy Scissors

Catalog Number	Tip	Tip Size (mm)	Blade	Handle	Length (in.)
EL-4883	B/B	10	ANG, CVD	FLAT, WIDE SERR	4 3/4

GILLS-WELSH VANNAS Scissors


Catalog Number	Tip	Tip Size (mm)	Handle	Length (in.)
EL-4884	S/S	11	FLAT, WIDE SERR	3 1/4

DEWECKER Iris Scissors


Catalog Number	Tip	Blade	Length (in.)
EL-4885	S/B	CVD	4 1/2



Eye Scissors



EL-4886




EL-4886



EL-4887-EL-4888




EL-4887



EL-4888




EL-4889



EL-4889




EL-4890




EL-4890



EL-4891-EL-4892



EL-4891



EL-4892

WESTCOTT Stitch and Utility Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4886	S/S	CVD	FLAT, SERR	5 1/2
EL-4887	S/S	CVD	FLAT, SERR	4 1/2
EL-4888	B/B	CVD	FLAT, SERR	4 1/2

WESTCOTT Tenotomy Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4889	B/B	CVD, RIGHT HAND	FLAT, SERR	4 3/4
EL-4890	B/B	CVD, LEFT HAND	FLAT, SERR	4 3/4

MILLESI Micro Scissors

Catalog Number	Tip	Blade	Handle	Length (in.)
EL-4891	S/S	STR SERR	FLAT, SERR	6 1/4
EL-4892	S/S	CVD, SERR	FLAT, SERR	6 1/4




EL-4896 $\frac{1}{2}$




EL-4896 $\frac{1}{1}$



EL-4895-EL-4897 $\frac{1}{2}$




EL-4897 $\frac{1}{1}$




EL-4895 $\frac{1}{1}$



EL-4893-EL-4898 $\frac{1}{2}$




EL-4893 $\frac{1}{1}$



EL-4898 $\frac{1}{1}$




EL-4894 $\frac{1}{2}$



EL-4894 $\frac{1}{1}$




EL-4900 $\frac{1}{2}$




EL-4900 $\frac{1}{1}$



EL-4901-EL-4903 $\frac{1}{2}$




EL-4901 $\frac{1}{1}$



EL-4903 $\frac{1}{1}$



EL-4902 $\frac{1}{2}$



EL-4902 $\frac{1}{1}$



EL-4904 $\frac{1}{2}$



EL-4904 $\frac{1}{1}$

BARRAQUER Micro Needle Holder, Round Knurled Handle W/O Catch













Catalog Number	Tip	Length (in.)
EL-4893	STR	4 3/4
EL-4894	STR	5 1/4
EL-4895	CVD, DEL	4 3/4
EL-4896	CVD	4 1/4
EL-4897	CVD	4 1/2
EL-4898	CVD	4 3/4
EL-4900	CVD	5 1/4

BARRAQUER Micro Needle Holder, Round Knurled Handle W/ Catch

Catalog Number	Tip	Length (in.)
EL-4901	STR	4 3/4
EL-4902	STR	5 1/4
EL-4903	CVD	4 3/4
EL-4904	CVD	5 1/4



Eye Micro Needle Holders



COHEN Micro Needle Holder, Round Knurled Handle

Catalog Number	Tip	Handle	Length (in.)
EL-4905	CVD	W/O CATCH	4 1/8
EL-4906	CVD	W/ CATCH	4 1/8

TROUTMAN Micro Needle Holder, Flat Handle W/ Catch

Catalog Number	Tip	Length (in.)
EL-4907	CVD	4 3/4
EL-4908	CVD	5 1/8

CASTROVIEJO Micro Needle Holder, Flat Handle

Catalog Number	Tip	Jaw Length (mm)	Handle	Length (in.)
EL-4909	STR	9	W/ CATCH	5 1/2
EL-4910	STR	9	W/ CATCH	5 1/8
EL-4911	STR	16	W/ CATCH	5 7/8
EL-4912	STR	10	W/ CATCH	5 3/8
EL-4913	STR	10	W/O CATCH	5 3/8
EL-4914	CVD	10	W/O CATCH	5 3/8


Micro Needle Holder, Round Knurled Handle W/O Catch

Catalog Number	Tip	Length (in.)
EL-4915	DEL, CVD	4 3/4

Eye Micro Needle Holders



EL-4916-EL-4917



EL-4916



EL-4917

Micro Needle Holder, Round Knurled Handle W/ Catch

Catalog Number	Tip	Length (in.)
EL-4916	STR	6
EL-4917	CVD	6



EL-4918




EL-4918

HAMPTON Micro Needle Holder, Heavy Jaw, Round Knurled Handle W/ Catch


Catalog Number	Tip	Length (in.)
EL-4918	CVD	5 1/8



EL-4919-EL-4921




EL-4919




EL-4921



EL-4920-EL-4922



EL-4920




EL-4922

MCPHERSON Micro Needle Holder, Flat Handle, Delicate


Catalog Number	Tip	Handle	Length (in.)
EL-4919	STR	SERR W/ CATCH	4 1/4
EL-4920	STR	FENESTRATED W/O CATCH	4 1/4
EL-4921	CVD	SERR W/O CATCH	4 1/8
EL-4922	CVD	FENESTRATED W/ CATCH	4 1/8



Eye Forceps




EL-4923



EL-4923




EL-4924



EL-4924



EL-4925



EL-4925



EL-4926



EL-4926




EL-4927



EL-4927



EL-4928



EL-4928



EL-4929




EL-4929

Micro Jewelers Forceps

Catalog Number	Tip	Tip Size (mm)	Tip Size (#)	Length (in.)
EL-4923	STR, W/ PLATFORM, X-SERR	3	-	4 1/4
EL-4924	STR, FINE, W/ PLATFORM	6	1	4 3/4
EL-4925	STR, FINE	-	1	4 3/4
EL-4926	STR, X-FINE, W/ PLATFORM	5	3	4 3/8
EL-4927	STR	0.3	3	5 1/4
EL-4928	STR, X-FINE, W/ PLATFORM	5	3	4 3/4
EL-4929	STR, X-FINE, W/ PLATFORM	6	4	4 3/8
EL-4930	STR	0.2	4	4 1/4
EL-4931	STR, X-FINE, W/ PLATFORM	6	5	4 3/8
EL-4932	STR	0.15	5	4 1/4
EL-4933	CVD	0.20	5	4 1/2




EL-4930



EL-4930




EL-4931



EL-4931



EL-4932



EL-4932



EL-4933



EL-4933




EL-4934



EL-4934



EL-4935



EL-4935



EL-4936



EL-4936



EL-4937



EL-4937




EL-4938



EL-4938



EL-4939-EL-4941



EL-4939




EL-4940



EL-4941




EL-4942-EL-4944



EL-4942



EL-4943



EL-4944

Micro Forceps

Catalog Number	Tip	Jaw	Length (in.)
EL-4934	STR	SMOOTH	4
EL-4935	ANG	SMOOTH	4

BARRAQUER Cilia Forceps

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-4936	SLIGHT-ANG	6	4 1/9

WISE Micro Suture Forceps W/ Platform

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-4937	STR	0.3	5 1/4
EL-4938	CVD	0.3	5 1/4





Iris Forceps, Delicate

Catalog Number	Tip	Jaw	Length (in.)
EL-4939	STR	SERR	4
EL-4940	CVD	SERR	4
EL-4941	HEAVY-CVD	SERR	4

Iris Forceps, X-Fine

Catalog Number	Tip	Teeth	Length (in.)
EL-4942	STR	1 x 2	4
EL-4943	CVD	1 x 2	4
EL-4944	HEAVY-CVD	1 x 2	4

Eye Forceps



BISHOP-HARMON Iris Forceps

Catalog Number	Teeth	Jaw	Length (in.)
EL-4945	1 x 2	-	3 1/4
EL-4946	-	SERR	3 1/2

ST. MARTIN Tissue Forceps, Delicate

Catalog Number	Teeth	Length (in.)
EL-4947	1 x 2	3 1/8

CASTROVIEJO Tissue Forceps, Delicate

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4948	0.9	1 x 2	4


CASTROVIEJO Fixation Forceps, W/O Platform

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4949	0.5	1 x 2	4 1/4




**CASTROVIEJO Suture Forceps,
W/ 1 x 2 Teeth**


Catalog Number	Tip	Tip Size (mm)	Platform	Handle	Length (in.)
EL-4950	STR	0.12	YES	SERR	4 1/8
EL-4951	STR	0.12	YES	SERR	4 3/8
EL-4952	STR	0.22	YES	SERR	4 1/4
EL-4953	STR	0.3	YES	SERR	4 1/8
EL-4954	STR	0.3	YES	SERR	4 3/8
EL-4955	STR	0.5	YES	WIDE	4 1/4
EL-4956	STR	0.5	YES	SERR	4 1/4
EL-4957	STR	0.5	NO	WIDE W/ LOCK	4 1/4
EL-4958	STR	0.9	YES	SERR	4 1/4
EL-4959	STR	1.5	YES	WIDE	4 1/4
EL-4960	ANG 45°	0.5	YES	SERR	4 1/4



Eye Forceps



EL-4961



EL-4961




EL-4962



EL-4962




EL-4963



EL-4963



EL-4964



EL-4964



EL-4965-EL-4966



EL-4965



EL-4966



EL-4967



EL-4967

CASTROVIEJO Capsule Forceps, Cross-Action

Catalog Number	Tip	Length (in.)
EL-4961	10 MM FROM ANGLE TO TIP	4 1/8

CASTROVIEJO Colibri Type Corneal Forceps

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4962	0.12	1 x 2	4 1/8

BONN Eye Forceps, W/ Platform

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4963	0.2	1 x 2	3 3/4

GRAEFE Tissue Forceps




Catalog Number	Length (in.)
EL-4964	4 1/4

Suture Tying Forceps, W/ Platform


Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-4965	STR	0.3	4 3/4
EL-4966	CVD	0.3	4 3/4

OSHEL-ELSCHNIG Forceps, W/ Platform

Catalog Number	Tip	Tip Size (mm)	Teeth	Length (in.)
EL-4967	0.12	3	1 x 2	4 1/8



EL-4971



BONACCOLTO Utility Forceps

Catalog Number	Tip Size (mm)	Length (in.)
EL-4968	1.2 x 16	4 1/4
EL-4969	1.7 x 17	4 1/4

PUTTERMAN-Type Muscle Forceps, W/ Lock

Catalog Number	Tip	Jaw	Length (in.)
EL-4970	6 PINS	SERR	4

TENNANT Tying Forceps, W/ Platform


Catalog Number	Tip	Tip Size (mm)	Jaw	Handle	Length (in.)
EL-4971	X-DEL	6	SMOOTH	ROUND	4 1/4

SINSKEY Tying Forceps

Catalog Number	Tip	Jaw	Handle	Length (in.)
EL-4972	X-DEL	SMOOTH, CVD	FLAT, WIDE SERR	4




Eye Forceps



EL-4973

1/2



EL-4973

1/1




4/1



EL-4974

1/2



EL-4974

1/1




2/1



EL-4975

1/2



EL-4975

1/1



1/1




1/1



EL-4976

1/2



EL-4976

1/1




4/1



EL-4977

1/2



EL-4977

3/1

1/1

SHEPARD Lens Holding Forceps, Cross-Action

Catalog Number	Tip	Jaw	Handle	Length (in.)
EL-4973	DUCKBILL	CVD	FLAT, WIDE SERR	5

CLAYMAN Lens Holding Forceps, Cross-Action

Catalog Number	Tip	Jaw	Handle	Length (in.)
EL-4974	DEL, ANG	CVD	FLAT, WIDE SERR	4 3/4

Scleral Plug Forceps, Cross-Action

Catalog Number	Tip	Tip Size (mm)	Jaw	Handle	Length (in.)
EL-4975	CVD	9	ANG	FLAT, SERR	4 1/8

DEVGAN Titanium Capsulorhexis Forceps, Cross-Action

Catalog Number	Tip	Handle	Length (in.)
EL-4976	X-DEL, SEMI-SHARP	ROUND	5


GRIFFITHS-BROWN Tissue Forceps

Catalog Number	Tip	Length (in.)
EL-4977	DEL	4 1/4



EL-4978

1/1




EL-4978



EL-4979

1/1



EL-4979



EL-4980

1/1



EL-4980



EL-4981

1/1



EL-4981

4/1




EL-4982-EL-4983

1/1



EL-4982



EL-4983

BONN Suture Forceps, W/ Platform

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4978	0.12	1 x 2	3 3/4

BONN Tissue Forceps

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4979	0.12	1 x 2	3 3/4

THORPE Conjunctival Fixation Forceps, W/ Platform

Catalog Number	Teeth	Length (in.)
EL-4980	2 x 3	3 3/4

THORPE Corneal Forceps W/ Platform




Catalog Number	Tip Size (mm)	Jaw	Teeth	Length (in.)
EL-4981	0.12	STR	2 x 3, ANG 45°	3 3/4

LESTER Fixation Forceps

Catalog Number	Teeth	Length (in.)
EL-4982	FINE, 1 x 2	3 3/4
EL-4983	FINE, 2 x 3	3 3/4



Eye Forceps



NUGENT Utility Forceps

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
EL-4984	SERR, ANG 45°	10	4 1/4
EL-4985	SMOOTH, ANG 45°	10	4 1/4

TROUTMAN Superior Rectus Forceps

Catalog Number	Tip Size (mm)	Jaw	Teeth	Length (in.)
EL-4986	10	ANG	1 x 2	4 1/8

BERKE Ptosis Forceps, W/ Lock


Catalog Number	Jaw Length (mm)	Length (in.)
EL-4987	20	4
EL-4988	27	4 1/8

Serreŉne Clamp

Catalog Number	Tip	Jaw	Jaw Length (mm)	Length (in.)
EL-4989	STR	SERR	10	1 3/8

LEHNEL-UTRATA Capsulotomy Forceps

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-4990	ANG	12	3 1/2



UTRATA Capsulorhexis Forceps

Catalog Number	Tip	Tip Size (mm)	Handle	Length (in.)
EL-4991	ANG	12	FLAT, SERR	3 1/8
EL-4992	ANG	11	FLAT, FENESTRATED	4

COLIBRI-Style Corneal Forceps, W/ Platform

Catalog Number	Tip	Tip Size (mm)	Teeth	Length (in.)
EL-4993	ANG 50°	6	1 x 2	4 1/4

TROUTMAN-BARRAQUER Corneal Forceps

Catalog Number	Tip Size (mm)	Teeth	Length (in.)
EL-4994	0.12	1 x 2	2 7/8

Titanium COLIBRI Forceps



Catalog Number	Tip Size (mm)	Jaw	Teeth	Length (in.)
EL-4995	0.12	SMOOTH	1 x 2	3

MCPHERSON Angled Corneal Suturing Forceps, W/ Platform

Catalog Number	Tip	Tip Size (mm)	Teeth	Length (in.)
EL-4996	ANG 45°	5	1 x 2	3 3/8



Eye Forceps



BRACKEN Iris Forceps

Catalog Number	Tip	Jaw	Length (in.)
EL-4997	CVD	CROSS-SERR	3 7/8

KELMAN-MCPHERSON Forceps, W/ Platform

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-4998	ANG 45°	8.5	3 3/8

MCPHERSON Forceps, W/ Platform



Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-4999	STR	5	3 1/4

BECHERT-MCPHERSON Forceps

Catalog Number	Tip	Tip Size (mm)	Jaw	Length (in.)
EL-5000	ANG	10	SMOOTH	3 1/4

SHEETS-MCPHERSON Forceps

Catalog Number	Tip	Tip Size (mm)	Jaw	Length (in.)
EL-5001	1 MM X-SERR	10	ANG	3 1/4



JAMESON Muscle Forceps, W/ Lock

Catalog Number	Tip Size (mm)	Jaw	Length (in.)
EL-5002	12	ANG, LEFT	3 3/8
EL-5003	12	ANG, RIGHT	4

MCCULLOUGH Utility Forceps

Catalog Number	Tip	Tip Size (mm)	Jaw	Length (in.)
EL-5004	STR	15	SMOOTH	3 5/8

Tying Forceps


Catalog Number	Jaw	Handle	Length (in.)
EL-5005	SMOOTH	FLAT, SERR	4 1/4

Alabama University Forceps, W/ Platform

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-5006	V-GROOVED	15	3 1/2



Eye Forceps



Von GRAEFE Iris Forceps

Catalog Number	Tip	Teeth	Length (in.)
EL-5007	STR	1 x 2	3

GILL Iris Forceps

Catalog Number	Tip	Jaw	Length (in.)
EL-5008	CVD	CROSS-SERR	3


WILLS Hospital Utility Forceps

Catalog Number	Tip	Jaw	Length (in.)
EL-5009	STR	SERR	4

WATZKE Sleeve Spreading Forceps

Catalog Number	Tip	Jaw Length (mm)	Length (in.)
EL-5010	ANG	10	4 3/4







EL-5012





EL-5013




EL-5014



EL-5011



EL-5015



EL-5016

BARRAQUER Eye Speculum

Catalog Number	Blade Size (mm)	Frame	Length (in.)
EL-5011	10	WIRE	1 1/2
EL-5012	10	WIRE	1 3/4
EL-5013	14	HEAVY WIRE	1 3/4
EL-5014	15	SOLID	1 3/4

KRATZ-BARRAQUER Eye Speculum

Catalog Number	Blade Size (mm)	Blade Shape	Length (in.)
EL-5015	10	PEDIATRIC	1 3/8
EL-5016	14	ADULT	1 1/2



Eye Specula



EL-5017



EL-5018



EL-5019



EL-5020



EL-5021



EL-5022



EL-5023



EL-5024

ALPHONSO Speculum, Infant


Catalog Number	Blade Size (mm)	Length (in.)
EL-5017	7 x 5 AND 9 WIDE	1 3/4

LIEBERMAN Speculum, Adjustable

Catalog Number	Blade Size (mm)	Blade Shape	Length (in.)
EL-5018	16	STR, SQUARE WIRE	3 1/8
EL-5019	14	ANG, SQUARE WIRE	3 1/4
EL-5020	14	ANG, V-SHAPE WIRE	3 1/4
EL-5021	14	ANG, SOLID	3 1/4
EL-5022	14	ANG, REVERSIBLE, SOLID	3 1/8


JAFFE Wire Lid Retractor, Pair

Catalog Number	Blade Size (mm)	Length (in.)
EL-5023	13	1 1/2
EL-5024	15	1 1/2





EL-5025



EL-5025

SAUER Eye Speculum

Catalog Number	Length (in.)
EL-5025	1 3/8



EL-5026




EL-5026

LANCASTER Eye Speculum


Catalog Number	Length (in.)
EL-5026	2 3/4
EL-5027	3 1/8




EL-5027



EL-5027



EL-5028




EL-5028

LOGAN Modified REILL Self-Retaining Retractor

Catalog Number	Length (in.)
EL-5028	2



Eye Hooks and Retractors



EL-5029-EL-5032



EL-5029



EL-5031

SEWALL Orbital Retractor

Catalog Number	Blade Size (mm)	Length (in.)
EL-5029	7 x 52	6 3/4
EL-5030	9 x 44	6 3/8
EL-5031	11 x 50	6 3/8
EL-5032	13 x 67	6 3/4



EL-5030



EL-5032



EL-5033-EL-5037

DESMARRES Lid Retractor

Catalog Number	Blade Size (mm)	Length (in.)
EL-5033	9 x 9	5 1/2
EL-5034	9 x 12	5 1/2
EL-5035	11 x 14	5 1/2
EL-5036	12 x 16	5 1/2
EL-5037	13 x 18	5 1/2



EL-5038

SCHEPENS Orbital Retractor

Catalog Number	Tip	Blade Size (mm)	Length (in.)
EL-5038	SMOOTH	14 x 61	6 1/8



EL-5039


JAEGER Lid Plate, Clear Plastic

Catalog Number	Blade Size (mm)	Length (in.)
EL-5039	20/23	3 3/4



Eye

Hooks / Retractors / Manipulators



BERKE-JAEGER Lid Plate

Catalog Number	Blade Size (mm)	Length (in.)
EL-5040	20/30	4 $\frac{1}{4}$

WELLS Enucleation Spoon

Catalog Number	Diameter (mm)	Length (in.)
EL-5041	21	5 $\frac{3}{4}$

BUNGE Evisceration Spoon

Catalog Number	Tip	Diameter (mm)	Length (in.)
EL-5042	SMALL	8	5 $\frac{1}{4}$
EL-5043	LARGE	12	5 $\frac{1}{2}$

SCHOCKET Scleral Depressor, Double Ended

Catalog Number	Length (in.)
EL-5044	5 $\frac{1}{4}$


FLYNN Scleral Depressor

Catalog Number	Diameter (mm)	Length (in.)
EL-5045	5	3 $\frac{1}{8}$




Eye

Hooks / Retractors / Manipulators




WHITMAN-SHEPHERD Capsule Polisher, Double Ended

Catalog Number	Tip	Length (in.)
EL-5046	LEFT ANG, RIGHTANG	5 1/8




KOCH Stop and Chop Manipulator, Angled Shaft

Catalog Number	Tip Size (mm)	Length (in.)
EL-5047	1.5	4 3/4



CAMELLIN Style LASIK Alcohol Well

Catalog Number	Diameter (mm)	Length (in.)
EL-5048	8.5	4 1/2





SLOANE LASIK Spatula

Catalog Number	Tip	Description	Length (in.)
EL-5049	SEMI-CIRCULAR, BEVELED EDGE	MICRO HOE	4 1/2
EL-5050	SEMI-SHARP	EPI PEELER	4 1/2



Eye

Hooks / Retractors / Manipulators



CHANG-SEIBEL Phaco Chopper, Double Ended

Catalog Number	Tip	Style	Length (in.)
EL-5051	ANG C-SHAPED BLUNT CHOPPER / VERTICAL SAFETY CHOPPER	F/ LEFT HAND	5 1/2

KOCH Nucleus Spatula

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-5052	SLIGHTLY CVD W/ DUCKBILL TIP	1	4 1/2

LIEBERMAN Nucleus Manipulator

Catalog Number	Style	Shaft	Length (in.)
EL-5053	F/ LEFT HAND	ANG	4 1/2

CHANG Combo Phaco Chopper, Double Ended


Catalog Number	Tip	Style	Shaft	Length (in.)
EL-5054	C-SHAPED BLUNT CHOPPER / 1.3 MM SHARP QUICK CHOP TIP	F/ LEFT HAND	ANG	5 1/4

NAGAHARA Nucleus Spatula

Catalog Number	Tip	Style	Shaft	Length (in.)
EL-5055	FLAT INTERIOR EDGE	F/ LEFT HAND	ANG	4 1/4

Eye

Hooks / Retractors / Manipulators



BECHERT Nucleus Rotator

Catalog Number	Tip	Shaft	Length (in.)
EL-5056	Y-SHAPED	ANG 45° 7 MM FROM TIP	4 3/4

JAFFE-BECHERT Nucleus Rotator

Catalog Number	Tip	Shaft	Length (in.)
EL-5057	Y-SHAPED	ANG 45° 9 MM FROM TIP	4 3/4

DRYSDALE Nucleus Manipulator

Catalog Number	Tip	Length (in.)
EL-5058	PADDLE-SHAPED	4 1/4

PISACANO Nucleus Rotator

Catalog Number	Tip	Length (in.)
EL-5059	VERTICAL FLAT	4 1/4







ANIS-WAND Nucleus Manipulator

Catalog Number	Tip	Tip Size (mm)	Shaft	Length (in.)
EL-5060	BALL-SHAPED	0.5	STR	4 3/4



Eye

Hooks / Retractors / Manipulators



SIEBEL Flap Lifter and Spatula, Double Ended

Catalog Number	Tip	Length (in.)
EL-5061	Z-SHAPED AND 1.5 MM ANG POINT	4 1/8

JAFFE Lens Spatula

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-5062	FLAT	0.5	4 3/4

LESTER IOL Manipulator

Catalog Number	Tip	Tip Size (mm)	Shaft	Length (in.)
EL-5063	TEARDROP-SHAPED	0.25	STR	4 1/2
EL-5064	TEARDROP-SHAPED	0.25	ANG	4 1/2

LESTER II Titanium IOL Manipulator




Catalog Number	Tip	Tip Size (mm)	Shaft	Length (in.)
EL-5065	TEARDROP-SHAPED	0.20	ANG	4 1/2

MALONEY IOL Manipulator

Catalog Number	Tip	Tip Size (mm)	Shaft	Length (in.)
EL-5066	TEARDROP-SHAPED	0.25	CVD	4 1/2

Eye

Hooks / Retractors / Manipulators



BARRAQUER Cyclodialysis Spatula

Catalog Number	Blade Size (mm)	Blade	Length (in.)
EL-5067	0.5 x 15	ANG	4 $\frac{7}{8}$

CASTROVIEJO Cyclodialysis Spatula, Double Ended

Catalog Number	Blade Size (mm)	Length (in.)
EL-5068	0.5 x 10 / 0.5 x 15	5 $\frac{1}{8}$

CASTROVIEJO Synchia Spatula, Double Ended

Catalog Number	Blade Size (mm)	Length (in.)
EL-5069	0.5	5 $\frac{1}{4}$

NEW ORLEANS Lens Loop

Catalog Number	Tip	Length (in.)
EL-5070	4 x 10	5 $\frac{1}{2}$



PATON Spatula and Spoon, Double Ended

Catalog Number	Tip	Length (in.)
EL-5071	6 MM WIDE GROOVED SPOON / 1.5 MM FLAT BLADE ANG AT 8 MM	5 $\frac{1}{2}$



Eye

Hooks / Retractors / Manipulators



GRAETHER Collar Button

Catalog Number	Tip	Shaft	Length (in.)
EL-5072	0.6 MM BUTTON	ANG 45° 10 MM FROM TIP	4 3/4

CLAYMAN-JAFFE Lens Hook


Catalog Number	Tip	Tip Size (mm)	Shaft	Length (in.)
EL-5073	ANG 90° BLUNT	0.15	ANG 45° 10 MM FROM TIP	4 1/2

HARRS-SINSKEY Lens Manipulating Hook

Catalog Number	Tip	Tip Size (mm)	Shaft	Handle	Length (in.)
EL-5074	BLUNT	0.25	ANG 45° 10 MM FROM TIP	FLAT	4 1/2

SINSKEY Lens Manipulating Hook

Catalog Number	Tip	Tip Size (mm)	Shaft	Handle	Length (in.)
EL-5075	BLUNT	0.25	STR	FLAT	4 5/8
EL-5076	BLUNT	0.20	MODIFIED CVD	FLAT	4 1/2
EL-5077	BLUNT, REVERSE	0.20	ANG	FLAT	4 1/2
EL-5078	BLUNT	0.25	STR	ROUND	4 1/2
EL-5079	BLUNT	0.25	ANG	ROUND	4 1/2



SINSKEY II Lens Manipulating Hook

Catalog Number	Tip	Tip Size (mm)	Shaft	Handle	Length (in.)
EL-5080	BLUNT	0.2	STR	FLAT	4 1/2
EL-5081	BLUNT	0.2	ANG	FLAT	4 1/2
EL-5082	BLUNT	0.2	STR	ROUND	4 1/2
EL-5083	BLUNT	0.2	ANG	ROUND	4 1/2


SINSKEY II Titanium Lens Manipulating Hook

Catalog Number	Tip	Tip Size (mm)	Shaft	Handle	Length (in.)
EL-5084	BLUNT	0.2	STR	FLAT	4 1/2
EL-5085	BLUNT	0.2	ANG	ROUND	4 1/2



Eye

Hooks / Retractors / Manipulators



KUGLEN Iris Hook and Lens Manipulator

Catalog Number	Tip	Tip Size (mm)	Shaft	Handle	Length (in.)
EL-5086	K-HOOK	0.5	STR	SQUARE	4 5/8
EL-5087	PUSH-PULL	-	STR	FLAT	4 1/2
EL-5088	PUSH-PULL	-	ANG	FLAT	4 1/2
EL-5089	PUSH-PULL	-	STR	ROUND	4 1/2
EL-5090	PUSH-PULL	-	ANG	ROUND	4 1/2

KUGLEN Titanium Iris Hook and Lens Manipulator

Catalog Number	Tip	Shaft	Handle	Length (in.)
EL-5091	K-HOOK	STR	ROUND	4 1/2
EL-5092	PUSH-PULL	ANG	ROUND	4 1/2

Eye

Hooks / Retractors / Manipulators



EL-5093

1/2




EL-5093

4/1



EL-5094

1/2



EL-5094

4/1



EL-5095

1/2




EL-5095

4/1



EL-5096

1/2



EL-5096

2/1



EL-5097

1/2



EL-5097

4/1

OSHER "Y"-Hook Lens Manipulator

Catalog Number	Tip	Shaft	Length (in.)
EL-5093	Y-SHAPED	STR	4 1/2

JAFFE-MALTZMAN Lens Manipulating Hook

Catalog Number	Tip	Tip Size (mm)	Shaft	Length (in.)
EL-5094	V-SHAPED	0.15	ANG	4 1/2

MCPHERSON Iris Spatula

Catalog Number	Blade Size (mm)	Length (in.)
EL-5095	1 x 16	4 1/2

KIRBY Hook and Muscle Loop, Double Ended

Catalog Number	Blade Size (mm)	Length (in.)
EL-5096	5 x 10	6 1/8


BONN Micro Iris Hook

Catalog Number	Tip	Shaft	Length (in.)
EL-5097	BLUNT	STR	5



Eye

Hooks / Retractors / Manipulators



JAMESON Muscle Hook

Catalog Number	Size	Tip Size (mm)	Length (in.)
EL-5098	SMALL	9	5 1/4
EL-5099	LARGE	10	5 1/4

GREEN Strabismus Hook


Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-5100	FLAT	8	5 1/8

Twist Fixation Hook

Catalog Number	Tip	Length (in.)
EL-5101	RIGHT, CLOCKWISE TWIST	4 1/4

Eye

Hooks / Retractors / Manipulators



GASS Retinal Detachment Hook

Catalog Number	Length (in.)
EL-5102	5 5/8

STEVENS Tenotomy Hook

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-5103	CVD, BLUNT	5	4 5/8

ALFONSO Nucleus Trisector

Catalog Number	Tip Size (mm)	Length (in.)
EL-5104	1.5	4 3/4

TYRELL Iris Hook

Catalog Number	Tip	Tip Size (mm)	Length (in.)
EL-5105	SHARP	2	4 3/4



Eye

Fixation Rings and Markers



EL-5107-EL-5117



EL-5118



EL-5118



EL-5119



EL-5119

FLIERINGA Scleral Fixation Ring


Catalog Number	Diameter (mm)	Description
EL-5106	12-22	SET OF 8 RINGS
EL-5107	12	-
EL-5108	13	-
EL-5109	14	-
EL-5110	15	-
EL-5111	16	-
EL-5112	17	-
EL-5113	18	-
EL-5114	19	-
EL-5115	20	-
EL-5116	21	-
EL-5117	22	-

OSHER-NEUMANN Corneal Marker


Catalog Number	Tip	Length (in.)
EL-5118	8 BLADES	4 1/8

THORNTON Swivel Fixation Ring

Catalog Number	Tip	Diameter (mm)	Length (in.)
EL-5119	BLUNT TEETH	13	4 1/2



Eye Fixation Rings and Markers



MENDEZ Degree Gauge, Calibrated 10° from 0-180°

Catalog Number	Inner Diameter (mm)	Length (in.)
EL-5120	12	4 1/4

BORES Axis Marker

Catalog Number	Length (in.)
EL-5121	4

HENDERSON Degree Gauge

Catalog Number	Tip	Length (in.)
EL-5122	BEVELED FACE, SMOOTH	4 1/8

HENDERSON Alignment Marker

Catalog Number	Tip	Length (in.)
EL-5123	BLUNT PIN MARKINGS AT 3, 6 AND 9 O'CLOCK	4


HENDERSON Toric IOL Marker

Catalog Number	Tip	Length (in.)
EL-5124	2.5 MM AXIS MARKS	4



Eye

Fixation Rings and Markers



EL-5125-EL-5128



EL-5125



EL-5126



EL-5127



EL-5128



EL-5087-EL-5088



EL-5129



EL-5130



EL-5131



EL-5132




EL-5133

BORES Optical Zone Marker

Catalog Number	Inner Diameter (mm)	Length (in.)
EL-5125	3.5	4
EL-5126	6	4 1/8
EL-5127	7	4 1/8
EL-5128	8	4 1/8


HOFFER Optical Zone Marker

Catalog Number	Inner Diameter (mm)	Length (in.)
EL-5129	4.25	4
EL-5130	5	4
EL-5131	5.5	4
EL-5132	6	4
EL-5133	8	4 1/8




WILLIAMS Lacrimal Probe, Double Ended

Catalog Number	Tip Size (#)	Length (in.)
EL-5134	0000/000	5 1/8
EL-5135	00/0	5 1/2
EL-5136	1/2	5 1/2
EL-5137	3/4	5 1/2
EL-5138	5/6	5 1/2
EL-5139	7/8	5 1/2



BOWMAN Lacrimal Probe, Double Ended

Catalog Number	Tip Size (#)	Tip Size (mm)	Description	Length (in.)
EL-5140	-	-	SET OF 6	-
EL-5141	0000/000	-	-	5 5/8
EL-5142	00/0	0.7 / 0.8	-	5 5/8
EL-5143	1/2	0.9 / 1.1	-	5 5/8
EL-5144	3/4	1.3 / 1.4	-	5 5/8
EL-5145	5/6	1.5 / 1.6	-	5 5/8
EL-5146	7/8	1.8 / 1.9	-	5 5/8




Pigtail Lacrimal Probe, Double Ended, W/ Suture Holes

Catalog Number	Tip	Diameter (mm)	Length (in.)
EL-5147	BLUNT	8	5 1/8
EL-5148	BLUNT	8	6 1/8



Eye Probes and Dilators



EL-5149 1/2

2/1 2/1
EL-5149

EL-5150 1/2

2/1
EL-5150

EL-5151 1/2

2/1
EL-5151

EL-5152 1/2

2/1
EL-5152

EL-5153 1/2

2/1
EL-5153

CASTROVIEJO Lacrimal Dilator, Double Ended

Catalog Number	Tip	Length (in.)
EL-5149	FINE TAPPER/ MEDIUM TAPER	4 3/4

RUEDEMANN Lacrimal Dilator


Catalog Number	Tip	Length (in.)
EL-5150	12 MM TAPER	2 3/4

WILDER Lacrimal Dilator

Catalog Number	Tip	Tip Size (mm)	Size	Length (in.)
EL-5151	LONG TAPER	23	1	4
EL-5152	MEDIUM TAPER	18	2	4
EL-5153	SHORT TAPER	13	3	4

**CATALANO Lacrimal Intubation Set,
Single Use**


Catalog Number	Description
EL-5154	TWO 0.7 MM PROBES ATTACHED TO 250 MM LONG SILICONE TUBING



EL-5154

**JACKSON Lacrimal Intubation Set,
Single Use**

Catalog Number	Description
EL-5155	TWO 0.9 MM PROBES ATTACHED TO 250 MM LONG SILICONE TUBING




EL-5155



Eye Lacrimal Intubation Sets


CRAWFORD Lacrimal Intubation Set, Single Use


Catalog Number	Description
EL-5156	TWO OLIVE-TIPPED 0.8 MM PROBES ATTACHED TO 250 MM LONG SILICONE TUBING



CRAWFORD Retrieving Hook

Catalog Number	Compatible w/	Length (in.)
EL-5157	MA922R	5 1/8





Lacrimal Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5158	REINFORCED, MALLEABLE	23	STR	1 3/8
EL-5159	REINFORCED, MALLEABLE	23	CVD	1 3/8
EL-5160	-	23	STR	1

SHAHINIAN Lacrimal Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5161	0.3 MM SIDE PORT	23	STR	3/8

RANDOLPH Cyclodialysis Cannula


Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5162	THIN	19	ANG 13 MM FROM TIP	1 3/8
EL-5163	DEL	19	ANG 13 MM FROM TIP	1 1/4

TROUTMAN Olive Tip Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5164	1	25	SLIGHT-CVD	1



Eye Cannulas



JENSEN Capsule Polisher, Roughened Tip, W/ Front Opening

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5165	OLIVE	25	CVD	1
EL-5166	OLIVE	27	CVD	1

KNOLLE Capsule Polisher, Roughened Tip

Catalog Number	Tip	Tip Size (Gauge)	Length (in.)
EL-5167	2 MM ANG	30	3/4

BISHOP-HARMON Anterior Chamber Irrigator, 19 GA, Angled

Catalog Number	Description	Length (in.)
EL-5168	COMPLETE SET	-
EL-5169	SILICONE BULB, 28 x 55 MM	-
EL-5170	PLASTIC ADAPTER	-
EL-5171	CANNULA	1

Hydrodissection Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5172	FLATTENED, FRONT OPENING	25	ANG	3/4



BRACKEN Irrigating Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5173	FLATTENED BEVELED, FRONT OPENING	19	CVD	1 1/4



SLADE LASIK Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5174	FLATTENED SPATULA, FRONT OPENING	26	ANG	1/2







BURATTO Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5175	FLATTENED, 1 FRONT AND 2 SIDE PORTS	25	ANG	1



Eye Cannulas







MCINTYRE Anterior Chamber Cannula

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5177	BLUNT	23	STR	3/4
EL-5178	BLUNT	26	STR	3/4
EL-5179	BLUNT	30	STR	3/4
EL-5180	BLUNT	23	ANG 12 MM FROM TIP	3/4
EL-5181	BLUNT	26	ANG 12 MM FROM TIP	3/4
EL-5182	BLUNT	30	ANG 12 MM FROM TIP	3/4

SIMCOE Cortex Extractor, Side Port

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5183	SMOOTH, 0.3 MM	23	CVD	1/2
EL-5184	ROUGHENED, 0.3 MM	23	CVD	1/2
EL-5185	SMOOTH, 0.3 MM	23	CVD	3/4
EL-5186	ROUGHENED, 0.3 MM	23	CVD	3/4



MORRISON Aspirating Cannula, Front Opening

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5187	SMOOTH, BLUNT	23	ANG 7 MM FROM TIP	1

MCINTYRE-BINKHORST Cannula, Front Opening, Hook Shape

Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5188	SMOOTH, BLUNT	26	STR	3/4
EL-5189	SMOOTH, BLUNT	26	RIGHT ANG	3/4
EL-5190	SMOOTH, BLUNT	26	LEFT ANG	3/4

CHANG Nucleus Hydrodissector Cannula






Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5191	FLAT, BEVELED OPENING	27	RIGHT ANG	1/2

ATKINSON Retrobulbar Needle

Catalog Number	Tip Size (Gauge)	Shaft	Length (in.)
EL-5192	23	STR	1 3/8



Eye Cannulas



KELMAN Cystotome

Catalog Number	Tip Size (Gauge)	Shaft	Length (in.)
EL-5193	22	STR	1 1/4

KRATZ Capsule Scrubber





Catalog Number	Tip	Tip Size (Gauge)	Shaft	Length (in.)
EL-5194	ROUGHENED	23	ANG	3/4

Back—ush Infusion Cannula

Catalog Number	Tip Size (Gauge)	Shaft	Length (in.)
EL-5195	23	STR	1 3/8

NICHAMIN Hydrodissection Cannula

Catalog Number	Tip Size (Gauge)	Shaft	Length (in.)
EL-5196	25	ANG	1



Air Injection Cannula

Catalog Number	Tip Size (mm)	Tip Size (Gauge)	Shaft	Length (in.)
EL-5197	5	27	ANG	3/4
EL-5198	5	30	ANG	3/4

WEST Lacrimal Cannula, Side Port


Catalog Number	Tip Size (mm)	Tip Size (Gauge)	Shaft	Length (in.)
EL-5199	0.3	23	STR	1 1/4

SIMCOE Irrigating Aspiring Unit

Catalog Number	Tip Size (Gauge)	Description
EL-5200	23	ASPIRATION THROUGH 0.3 MM TOP PORT, IRRIGATION THROUGH FRONT PORT, USES 225 MM TUBING



Eye Miscellaneous



BRAUNSTEIN Fixed Caliper, Double Ended

Catalog Number	Tip Size (mm)	Length (in.)
EL-5201	3.0 / 3.5	2 $\frac{3}{4}$
EL-5202	3.5 / 4.0	2 $\frac{3}{4}$




CASTROVIEJO Eye Caliper

Catalog Number	Tip	Description	Length (in.)
EL-5203	STR	0 - 15 MM RANGE	3 $\frac{1}{8}$
EL-5204	STR	0 - 40 MM RANGE	6 $\frac{3}{4}$
EL-5205	CVD	0 - 20 MM RANGE	3 $\frac{3}{8}$

AKAHOSHI Combo Prechopper, Cross-Action

Catalog Number	Jaw	Jaw Length (mm)	Length (in.)
EL-5206	STR	4.5	4 $\frac{1}{4}$





Corneal Dissector Knife

Catalog Number	Blade	Shaft	Length (in.)
EL-5207	CVD	ANG	4 1/2

TOOKE Corneal Knife

Catalog Number	Blade	Blade Size (mm)	Length (in.)
EL-5208	STR	3 x 18	4 3/4

KELLY-DESCEMET'S Membrane Punctor

Catalog Number	Diameter (mm)	Length (in.)
EL-5209	1	5 1/8

CARTER Sphere Introducer

Catalog Number	Length (in.)
EL-5210	5

MEYERHOEFFER Chalazion Curette

Catalog Number	Tip Size (#)	Cup (mm)	Length (in.)
EL-5211	0	1.5	5
EL-5212	1	1.7	5
EL-5213	2	2	5
EL-5214	3	2.5	5
EL-5215	4	3.5	5





All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse, we retain the rights to recall our catalogues and price lists and to take legal actions.